

Technical Specification Datasheet

Product Name: Magnetic Stirrer with Timer
Product Code: DMS-2000
Packing: 1 piece and extra magnetic needle
Dimensions: 230mm × 180mm × 120mm
Weight: 2.5kg

Technical Specifications:

- Maximum Stirring Volume - 2000 ml
- Stirring Speed Range - 0-1250 rpm
- Timer range - 0-999 minutes
- Power Requirement - 220V/50Hz or 110V/60Hz
- Ambient Temperature - Room Temperature to 60 °C
- Relative Moisture – 80%
- Microprocessor based functions
- 316 SS top
- Large Custom LCD

Techinstro assures the above specification and dimensions of the product as per datasheet. We ensure the quality of the product with moisture proof packing. This sheet is also available on the website, where a researcher or purchaser can also download it merely by clicking on the pdf link. The terms and condition for the product purchase mentioned in the quotation. For safe product storage, tips and procedure provided with the material. For more queries, contact our technical support.