

Technical Datasheets (TDS)

- **Product Name** – Kynar Bag
- **Product Series** – Kyn-TI

Technical Specification:-

- CAS Number – NA
- HS Code – 39239090
- Stock number – 10090
- MOC – PVDF, Silicon, and Polypropelene (PP)
- Film Thickness – 3 mil
- Capacity – 0.5 L to 1000 L
- Colour – Transparent
- Form Factor – Flexible Solid
- Maximum Operating Temperature – 160° F (71° C)
- Storage – @ room temperature
- Delivery Mode – Air, surface or sea mode
- Customisation – Only for bulk quantity
- Available pack size – As per requirement
- Packing – Light Weight moisture proof

Techinstro assures the above specification and dimensions of the product as per the datasheet. We ensure the quality of the product with moisture-proof packing. This sheet is also available on the website, where researcher or purchaser can also download it merely by clicking on the pdf link. The terms and condition for the product purchase mentioned in the quotation. For safe product storage, tips and procedure provided with the material.

For more queries, contact our technical support