

Technical Specification Datasheet (TDS)

Product Name: Titanium Dioxide Nanoparticles (TiO₂)


Product Series: TiO₂-NANO

Properties:

Specifications	Values
Purity	99.9%
Average Particle Size	30 - 80 nm
SSA	150 m ² /g
Molecular Weight	79.8658 g/mol
Molecular Formula	TiO ₂
Melting Point	1843 °C
Bulk Density	0.35 g / cm ³
CAS No.	13463 - 67 - 7
Physical Form	Powder
Morphology	Near Spherical
Color	White

SEM IMAGES

.....



Techinstro assures the above specification and dimensions of the product as per the datasheet. We ensure the quality of the product with moisture-proof packing. This sheet is also available on the website, where a researcher or purchaser can also download it merely by clicking on the pdf link. The terms and conditions for the product purchase are mentioned in the quotation. For safe product storage, tips and procedures are provided with the material. For more queries, contact our technical